

# Portable Raman Acoustic Tweezers for contact-less, container-less $\mu$ Plastics detection

Pietro G. Gucciardi

Pillar Environment- Spoke 4 - CNR

SAMOTHRACE 2<sup>nd</sup> Year:  
Experimental Prototypes Demo Showcase



Finanziato  
dall'Unione europea  
NextGenerationEU




Ministero  
dell'Università  
e della Ricerca



Italiadomani  
PIANO NAZIONALE  
DI RIPRESA E RESILIENZA

- **Technological gap** in the **detection of  $\mu$ Plastics and nanoplastics ( $< 1 \mu\text{m}$ )** in aquatic ecosystems and in food.
- End users: **environmental and food safety operators** seek **fast, portable and sensitive tools** to assess the  $\mu$ Ps presence through the food chain, **especially for the nanoplastics fraction.**
- **Issues with existing techniques** (Micro-FT-IR / Raman – Pyrolysis-GC-MS)
  - Analysis of digested organic tissues requires **lab-scale equipment**
  - **Time-consuming scanning. Bottom/Lateral side not analyzed**
  - **Contamination and background from the container and the detection of the smallest micro- and nano- fractions.**
- Market: the global  **$\mu$ Plastics detection market** size accounted for USD 4.80 billion in 2024, grew to USD 5.10 billion in 2025 and is estimated to hit around **USD 6.5 billion by 2029**

THE PROBLEM TO BE SOLVED



Finanziato  
dall'Unione europea  
NextGenerationEU




Ministero  
dell'Università  
e della Ricerca



Italiadomani  
PIANO NAZIONALE  
DI RIPRESA E RESILIENZA

- OUR SOLUTION**
- **Acoustic Trapping** with low cost ultrasound transducers (as in cars sensors)
    - **Analyze particles levitated in air.** Solid (mm to 200  $\mu\text{m}$ ). Liquid droplets
    - **Translate and Rotate** structures around their vertical axis without touching
  - **Combination with Raman** enables the **chemical analysis** of particles or molecules and nanostructures dispersed in liquids
  - Unique advantages of **Portable Raman-AT**:
    - Container-less analysis  $\rightarrow$  **No contamination or bkg signals**
    - Contact-less rotation  $\rightarrow$  **360° analysis** of the sample
    - **Self-concentration** of samples dispersed in liquid droplets
  - **Example:** full surface analysis of a PS particle stained with fluorescent ink with no need of manual intervention



Finanziato  
dall'Unione europea  
NextGenerationEU





Ministero  
dell'Università  
e della Ricerca



Italiadomani  
PIANO NAZIONALE  
DI RIPRESA E RESILIENZA

# Portable Raman Acoustic Tweezers

- **Starting point:** Commercial AT + Labscale spectrometer (TRL 3)
- Implementation of a **portable spectrometer and computer**
- Setup of a portable Raman-AT assembly for **out-of-lab operation**
- **Validation** on the detection of micro and nanoplastics in Lab  
Detection of **microplastics in mussels' digestate**
- Automatic loading system, Illumination, Imaging
- Optical/Acoustic enclosure for **outdoor operation**
- First validation on the detection of microplastics in liquid, outdoor (TRL ~ 5)



Finanziato dall'Unione europea  
NextGenerationEU



Ministero dell'Università e della Ricerca

## NEXT STEPS UNDER SAMOTHRACE ECOSYSTEM

- **Validate the definitive system** (visit the boot) for microplastics detection in indoor/outdoor operation
- Assess the sensitivity (**particles concentration and size limits**) on the detection of micro and nanoplastics **in watery systems and in digested organic tissues (real systems)**
- Demonstrate the system **operating it in food and environmental science institutes** (TRL 6)
- Explore a **broader range of applications** (e.g. detection of chemical pollutants and adsorbates traces in water, textile fibres in air, ...)



Finanziato  
dall'Unione europea  
NextGenerationEU



Ministero  
dell'Università  
e della Ricerca



Italiadomani  
PIANO NAZIONALE  
DI RIPRESA E RESILIENZA




www.samothrace.eu



Consiglio Nazionale delle Ricerche Istituto per i Processi Chimico-Fisici

# VISIT OUR DEMO AND MEET US AT BOOTH N. 54




Silva



Nino



Maria Grazia



Pietro



Onofrio



Alex



NextGenerationEU europea



Ministero dell'Università e della Ricerca



Italiadomani PIANO NAZIONALE DI RIPRESA E RESILIENZA